

Item Name : Seperatory Funnel Plastic

Item Code : 919.2008.01

Item Description :

This is made of Polypropylene, giving it excellent contact clarity and strength. This pear shaped, autoclavable seperating funnel has a leakproof threaded screw cap. The Polypropylene stopcock, fitted with a self lubricating PTFE plug.

Material :Polypropylene



Item Image :

Cat. No.	Particulars	Particulars
919.2008.01	100 ml	2 Pcs
919.2008.02	250 ml	2 Pcs
919.2008.03	500 ml	2 Pcs