


**Item Name :** Retord Stand Hammer toned Painted

**Item Code :** 1103.2008.01

**Item Description :**

Made of cast iron Hammer toned Painted with Straight road adjusted for mounting perpendicularly to bases.

**Item Image :**



Cat. No.	Base size mm	Rod Length mm	Rod dia
1103.2008.01	150 x 100	450	10
1103.2008.02	175 x 125	600	12.5
1103.2008.03	200 x 125	600	12.5
1103.2008.04	200 x 125	600	12
1103.2008.05	225 x 150	750	12
1103.2008.06	250 x 150	600	12
1103.2008.07	250 x 200	750	12
1103.2008.08	300 x 175	900	12
1103.2008.09	250 x 150	750	12

# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>