


Item Name : Microscope Slides, Half White

Item Code : 222.2008.06

Item Description :

Optical glass, hydrolytic class 3. Thickness approx 1mm, Size Approx. 76x26mm. We recommend using slides with ground edges to reduce the risk of injury. To protect against penetrating humidity the entire packing unit is also optionally available sealed in an aluminium bag. Pack quantity 2500 = 50 boxes of 50

Item Image :



Cat. No.	Description
222.2008.06	Ground edges
222.2008.07	Ground edges, twin frosted end
222.2008.08	Cut edges
222.2008.09	Cut edges, twin frosted end