

**Item Name :** VDRL Slide KLine Concavity Slides

**Item Code :** 223.2008.01

**Item Description :**

White (Clear). Optical glass, hydrolytic class 3. Rectangular ground edges. Thickness 4 mm, size in mm: 75 x 56. Concavities of 16 mm dia, depth 1.75 mm. Pack of 1.

**Item Image :**



Cat. No.	Description
223.2008.01	One Concavity
223.2008.02	Two Concavities
223.2008.03	Three Concavities
223.2008.04	Four Concavity
223.2008.05	Six Concavity
223.2008.06	Ten Concavity
223.2008.07	Twelve Concavity

A/Cantt-133 006. Haryana, India.

+0171 – 2698107

[info@glacierindia.com](mailto:info@glacierindia.com)

<https://www.glacierindia.com>