

**Item Name :** Retort Rings Made Of Mild Steel

**Item Code :** 1106.2008.01S

**Item Description :**

Die pressed (Superior Quality) Chrome plated, with boss head for holding separating funnels.  
Also Available in white powder coated paint

**Item Image :**



Cat. No.	Size
1106.2008.01S	50
1106.2008.02S	62.5
1106.2008.03S	75
1106.2008.04S	100
1106.2008.05S	25
1106.2008.06S	150
1106.2008.07S	200

A/Cantt-133 006. Haryana, India.

+0171 – 2698107

[info@glacierindia.com](mailto:info@glacierindia.com)

<https://www.glacierindia.com>