

Item Name : Microscope Slide Cover

Item Code : 224.2008.01

Item Description :

Pure White (clear) Boro. Glass, hydraulic class 1, with excellent chemical resistance.
Thickness No.1 (0.13 to 0.17mm) Pack Quantities: Square shape: 2000 = 10 boxes of 200. Rectangular Shape: 1000 = 10 boxes of 100.

Item Image :



Cat. No.	Description	Size (mm)
224.2008.01	Square Shape	18x18
224.2008.02	"	20x20
224.2008.03	"	22x22
224.2008.04	"	24x22
224.2008.05	Rectangular Shape	24x40
224.2008.06	"	24x50
224.2008.07	"	24x60

44, Vikas puri, P.O. Indraprastha, Jagadhri Road, @glacierindia.com
A/Cantt-133 006. Haryana, India. <https://www.glacierindia.com>