

**Item Name :** Retort Clamp Cross Pattern

**Item Code :** 1122.2008.01

**Item Description :**

Made of aluminium alloy, Cross pattern white powder coated with MS rod to grip 15mm to 25mm circles. For holding burettes. AS per DIN 12894

**Item Image :**



Cat. No.	Type	Max. Temp.
1122.2008.01	Cork Lined	60
1122.2008.02	PVC dipped	80
1122.2008.03	Silicon Sleeved	280
1122.2008.04	Cork Lined	60