

Item Name : Allihn Condensers Bulb Condenser

Item Code : 229.2008.05

Item Description :



Item Image :

Cat. No.	Socket	Cone	Effective Length (mm.)
229.2008.05	14/23	14/23	160
229.2008.06	19/26	19/26	160
229.2008.07	14/23	14/23	250
229.2008.08	19/26	19/26	250
229.2008.09	24/29	24/29	250
229.2008.10	24/29	24/29	300
229.2008.11	29/32	29/32	300
229.2008.12	24/29	24/29	400
229.2008.13	29/32	29/32	400

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>