

Item Name : Liquid Expansion Apparatus

Item Code : 236A.2008.01

Item Description :

Demonstrates the different thermal expansion of various liquids thus simultaneous comparison of the volumetric expansion of upto five different liquids. The apparatus is used in various laboratories to analyze difference in expansion of various liquids over the same temperature range. This apparatus comprises of 5 glass graduated capillary tubes of equal bore & bulbs at the bottom, mounted on a support frame with metal trough for immersion of bulbs in hot water in it.



Item Image :

Cat. No.

236A.2008.01