

**Item Name :** Universal Clamp

**Item Code :** 1137.2009.01

**Item Description :**

Aluminum die casted, powder coated, chrome plated with M.S. steel rod. For holding all kinds of flask necks upto 50 - 65 mm circles. As per DIN 12894.

**Item Image :**



| Cat. No.     | Type            | Max Temp. |
|--------------|-----------------|-----------|
| 1137.2009.01 | Cork Lined      | 60        |
| 1137.2009.02 | PVC dipped      | 80        |
| 1137.2009.03 | Silicon Sleeved | 280       |