

Item Name : Gas Preparation kit

Item Code : 249.2008.01

Item Description :

Suitable for gas preparation experiments for producing gases in small volumes, which can be collected over water, such as hydrogen, oxygen, carbon dioxide etc. The kit comprises of 5 gas collection tubes about 150 x 25mm with solid stopper, gas preparation tube of size 150 x 25mm with side arm and one holed stopper, thistle funnel, Zshaped gas delivery tube, tubing connector, sn as a chimney shelf and a metal test tube stand. All components packed in a plastic box with lid for storage, which also functions as tank while gas preparation.

Item Image :



Cat. No.

249.2008.01