

**Item Name :** Condenser Clamp Mild Steel

**Item Code :** 1136.2008.01

**Item Description :**

Die pressed, felt lined black/white powder coated with black painted rod.

**Item Image :**



Cat. No.	Type Max.	Temp.	Dia (mm)
1136.2008.01	Cork Lined	60	40
1136.2008.02	PVC dipped	80	40
1136.2008.03	Silicon Sleeved	280	40
1136.2008.01	Cork Lined	60	50
1136.2008.02	PVC dipped	80	50
1136.2008.03	Silicon Sleeved	280	50
1136.2008.01	Cork Lined	60	65
1136.2008.02	PVC dipped	80	65

# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>