

Item Name : Density Hydrometers

Item Code : 401.2008.01

Item Description :

For Analytical Laboratories Shot & Wax Poised, Adjusted for use at 15°C M-50° SP. Max Length 270 mm, Accuracy ± 0.0002

Item Image :



Cat. No.	Range
401.2008.01	0.650 - 0.700 g/ml
401.2008.02	0.700 - 0.750 g/ml
401.2008.03	0.750 - 0.800 g/ml
401.2008.04	0.800 - 0.850 g/ml
401.2008.05	0.850 - 0.900 g/ml
401.2008.06	0.900 - 0.950 g/ml
401.2008.07	0.950 - 1.000 g/ml

A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>