

Item Name : Specific Gravity Beaume Glass Hydrometers

Item Code : 402.2008.01

Item Description :

For general purpose shot & wax poised, adjusted for use at 20°C

Item Image :



Cat. No.	SP. GR. Range	°Be Range
402.2008.01	Heavy Liquid	1.000 - 2.000 0 - 70°Be
402.2008.02	Light Liquid	0.700 - 1.000 10 - 70°Be
402.2008.03	1.000 - 1.500	0 - 50°Be
402.2008.04	1.500 - 2.000	50 - 70°Be
402.2008.05	Universal Hydrometer	0.700 - 2.000 70 - 70°Be