

Item Name : Pipette Pump

Item Code : 805.2008.01

Item Description :

Presenting a fast release pipetting device for precise pipetting & zip quickemptying. The Knurled thumbwheel can be rotated to draw up or dispense liquids while the releaser may be pushed in for zip- quick emptying. The moulded chuck has a silicon rubber collet inside to hold various pipettes of standard sizes. It resists acids & alkalies and can be easily disassembled for thorough cleaning. Material : Acrylonitvile Butandiene Styrene

Item Image :



Cat. No.	Particulars	Packing
805.2008.01	2 ml	12 Pcs
805.2008.02	10ml	12 Pcs
805.2008.03	25 ml	12 Pcs