

Item Name : Tripod Stand Base White Powder Coated

Item Code : 1157.2008.01

Item Description :

White stand tripod cast iron base white powder coated with mild steel chrome plated rod

Item Image :



| Cat. No. | Base Dia mm. | Rod height mm. | Rod Dia mm. |
|--------------|--------------|----------------|-------------|
| 1157.2008.01 | 175 | 600 | 12 mm |
| 1157.2008.02 | 275 | 600 | 12 mm |
| 1157.2008.03 | 375 | 900 | 16 mm |