

Item Name : Wash Bottles

Item Code : 816.2008.01

Item Description :

Glacier Wash Bottles are made of Low Density Polyethylene which gives these bottles a translucent & unbreakable character. Solutions can be poured easily by squeezing the bottle lightly. The cap is fitted with a flexible polyethylene delivery tube which can be aimed wherever required. Material : Polythylene

Item Image :



Cat. No.	Particulars	Packing
816.2008.01	125 ml	12 Pcs
816.2008.02	250ml	12 Pcs
816.2008.03	500ml	6Pcs
816.2008.04	1000ml	6Pcs