

Item Name : Separatory Funnel Holder

Item Code : 824.2008.01

Item Description :

This non-corrosive Separatory Funnel Holder non-corrosive Separatory Funnel Holder, moulded in Polypropylene, has a front opening which provides clear view of the solution and allows unobstructed placement of separating funnels. Two tapered wedges are provided to ensure firm grip on rods having dia 1/2" or 3/8". Material : Polypropylene

Item Image :



| Cat. No. | Particulars | Packing |
|-------------|--------------------------|---------|
| 824.2008.01 | Separatory Funnel Holder | 12 Pcs |