

**Item Name :** Animal Cage

**Item Code :** 827.2008.01

**Item Description :**

Available in two different sizes, these Cages provide excellent shelter for different animal species in the laboratory. The metallic top provides proper ventilation and has provisions for feeding the species. Material : Polypropylene / Stainless Steel

**Item Image :**



Cat. No.	Particulars	Packing
827.2008.01	290 x 220 x 140 mm.	1 Pcs
827.2008.02	430 x 270 x 150 mm.	1 Pcs
827.2008.03	430 x 270 x 150 mm.	12 Pcs
827.2008.04	Water Bottle - 250ml	12 Pcs