

Cat No.

2045.2008.01



GLACIER INDIA
(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)



Item Name : Infrared Spectrometer

Item Code : 2045.2008.01

Item Description :

Spectral Range: 7800 to 350 cm⁻¹

- Resolution: Better than 0.85cm⁻¹
- Wave Number: $\pm 0.01\text{cm}^{-1}$ • Beam splitter: Ge coated Kbr
- Scanning Speed: 5-step adjustable for different applications
- Signal to noise ratio : better than 15,000:1 (RMS value, at 2100cm⁻¹, resolution: 4cm⁻¹, detector: DTGS, 1 minute data collection)
- Infrared Source: Air-cooled, high efficiency Reflex • Detector: DTGS
- Data system: Compatible computer • IR Library : 11 IR libraries include
- Software : FT-IR software contains all routines needed for basic spectrometer operations, including library search, quantitation and spec

Item Image :



44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>