

Item Name : Connectors T & Y New

Item Code : 829.2008.01

Item Description :

These Connectors add to existing range of connectors. The crest & valley structure provided at all the ends of the connectors give excellent leakproof grip to the tube.
 Material : Polypropylene

Item Image :



Cat. No.	Particulars	Size	Packing
829.2008.01	T Connector	6mm	36 Pcs
829.2008.02	Y Connector	6mm	36 Pcs
829.2008.03	T Connector	8mm	36 Pcs
829.2008.04	Y Connector	8mm	36 Pcs
829.2008.05	T Connector	10mm	36 Pcs
829.2008.06	Y Connector	10mm	36 Pcs
829.2008.01	Cross Connector	6mm	36 Pcs
829.2008.02	Cross Connector	8mm	36 Pcs
829.2008.03	Cross Connector	10mm	36 Pcs

829.2008.01	L Connector	6mm	36 Pcs
829.2008.02	L Connector	8mm	36 Pcs
829.2008.03	L Connector	10mm	36 Pcs

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>