

**Item Name :** Petri Dish Disposable, Sterilized

**Item Code :** 849.2008.01

**Item Description :**

These transparent, ready to use, disposable Petri Dishes are made of optically clear, nontoxic, Polystyrene, These Petri Dishes packed in sterilized by radiation. Material : Polystyrene

**Item Image :**



Cat. No.	Particulars	Packing
849.2008.01	90 mm. Radiation Stz.	460 Pcs
849.2008.02	60 mm. Radiation Stz.	300 Pcs
849.2008.03	35 mm. Radiation Stz.	500 Pcs