


Item Name : Cuvettes

Item Code : 850.2008.01

Item Description :

Acuvette is a kind of laboratory plasticware, usually a small tube of circular or square cross shape, closed at one end and open at the other, made of plastic made to hold samples for spectroscopic experiments. The best cuvettes are as clear as possible, without impurities that might affect a spectroscopic reading. Like a test tube, a cuvette may be open to the air on top or have a glass or other material to close it. Material : Polystrene

Item Image :



Cat. No.	Particulars	Packing
850.2008.01	Student Se	1 Pcs
850.2008.02	Teacher Set	1 Pcs