


Item Name : MCT-Twin Rack

Item Code : 860.2008.01

Item Description :

Provides this Twin Purpose Rack for both 0.5 ml & 1.5 ml. micro centrifuge tubes. Each side has 108 holes for 0.5 & 1.5 ml. MCT'S, respectively. The transparent top provides an excellent view as well as a protection cover for the tubes stored. Each face of the base has a unique alpha- numeric index for easy identification of the tube. Material : Polypropylene

Item Image :



Cat. No.	Particulars	Packing
860.2008.01	Rack for 0.5ml & 1.5 ml. MCTS	4Pcs