

**Item Name :** Constant Volume Air Thermometer

**Item Code :** 3046.2008.01

**Item Description :**

Consists of a polished teakwood stand approx. 750mm high fitted with a wooden scale reading from 0 to 60cm x 1mm & supporting an adjustable mercury reservoir by a metal rod. Adetachable air thermometer bulb about 70mm diameter connected via three way stopcock to mercury tube with zero ind surrounding the bulb supported by detachable metal shelf mounted on insulated asbestos screen. Mounted on a teakwood stand with a deep tray, which serves to catch.

**Item Image :**



Cat. No.

3046.2008.01