

**Item Name :** Physical Weight Boxes

**Item Code :** 1223.2008.01

**Item Description :**

Made Of Polished Brass ,Comes Nicely Packed In Velvet Lined Wooden Case With Hined Lid. Individual Masses Rest In Separate Compartments And Includes Small Covered Box Which Holds Nickel Silver And Aluminum Fractions, With Chrome Plated Brass Forceps To Lift Masses. Each Set Consists Of Masses In 1-2-2-5 System In The Desired Range.

**Item Image :**



Cat. No.

1223.2008.01