

Item Name : Characterstics Of PNP Transistors

Item Code : 3068.2008.01

Item Description :

Characteristics of a PNP : This is a complete unit comprising of two independent power supplies for base-emitter and collector-emitter junction connections of power supplies, meters and transistor with the help of patch cords .Rugged construction compact unit with self explanatory

Item Image :



Cat. No.

3068.2008.01

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>