

**Item Name :** Friction Board Apparatus

**Item Code :** 1344.2008.01

**Item Description :**

For experiments in friction and its dependence on various factors.

Comprises a wooden friction board about 490x75x18mm in size, a long aluminium friction board that fits over the wooden friction board, a wooden alider about 140x70x18mm with hook, a plain minium slider that fits over wooden slider and an aluminium slider with rubber sheet pasted on one side to provide various combinations of surfaces for friction experiments.

**Item Image :**



Cat. No.

1344.2008.01