

Item Name : Conservation of Linear Momentum

Item Code : 3069.2008.01

Item Description :

Versatile unit allowing the study of Conservation of Linear momentum

- To be a versatile apparatus which covers other kinetic principles such as motion, acceleration, Newton's laws of motion, velocity
- To comprise of two trolleys which runs on a precision air track which can be inclined to the horizontal
- Trolleys projected along air track using catapult mechanism
- Trolley mass to be adjusted by adding weights
- To be supplied with timing gates and timer
- Comprehensive technical manual • Set of slotted weights included



Item Image :

Cat. No.

3069.2008.01