

Item Name : Demonstration Meter Small Wooden

Item Code : 1346.2008.01

Item Description :

A useful meter for demonstration of the functioning of simple, permanent magnet type moving coil meters to a group of students. Movement with spring control and magnetic damping mounted inside a rectangular wooden frame with clear glass window at the front and ground glass and electrical connections are visible. Long pointer needle with thick arrowhead provides clear visibility from a distance. With 4mm socket terminals on one side for connecting input. Has zero adjustment screw at the front. Accuracy class index-2.5.

Item Image :



| Cat. No. | Size |
|--------------|-------|
| 1346.2008.01 | Big |
| 1346.2008.01 | Small |