

Cat No.

2064.2008.01



**GLACIER INDIA**  
**(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)**



**Item Name :** Hot Plate (Round)

**Item Code :** 2064.2008.01

**Item Description :**

Our hot plates are offered with three heat control switch, round plate is made of cast iron and body is made of thick mild steel. Rating 1.5 kw, works on 230 volts AC/DC.  
Size Available:- • 8" Dia • 12" Dia

**Item Image :**



# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>