

**Item Name :** Hot Plate With Power Cord

**Item Code :** 3078.2008.01

**Item Description :**

Laboratory hot plates are available in round and rectangular shapes.

Body is fabricated out of thick mild steel sheet duly finished in stoving enamel paint. The hot plate made of cast iron, nicely & Precisely machined & smoothed duly finished in heat resistant black paint is firmly mounted on the body. Heavy duty heating elements are securely layed under the plate to operate on 220 V AC 50 Hz single phase. Temp. is controlled by a three-heat rotary switch along with an indicator lamp for the mains and cord & plug.

A. HOT PLATE : (Round) Single - with 20 cm.(8") dia. Moulded top.

B. HOT PLATE : (Round) Double- with two plates of 20 cm (8") dia. having independent pilot lamps & switches

Optional:-

- Sunvic energy regulator can be provided against routine energy regulator.
- Digital temp. Indicator- cum- controller

**Item Image :**



Cat. No.

3078.2008.01