

**Item Name :** Hypsometer Copper

**Item Code :** 3079.2008.01

**Item Description :**

For determining the upper fixed point of the thermometers and boiling points at barometer pressure other than atmospheric. Comprising an upper boiling vessel with a steam jacketed inner tube. A lid with a tubule is provided which holds the thermometer under test in a bung. A side outlet draining tube and manometric tube mounted in a tubular with rubber bung. Overall Size:- 300 x 100 mm diameter.

Cat

**Item Image :**



Cat. No.

3079.2008.01