

**Item Name :** Inclined Plane With Brass Pulley

**Item Code :** 3081.2008.01

**Item Description :**

Inclined red wood plane 600 x 75 mm with pulley at one end, the other being hinged to wood base. A metal linear scale is fixed to edge of plane & a vertical linear and an angular 0-45° scale is fixed to base, with a clamp which also serves as an index. With a metal roller 60 x 25 mm length x dia., a wooden block and a Brass scale pan with cord.

**Item Image :**



Cat. No.

3081.2008.01