

**Item Name :** Slotted Weight Iron With Hanger

**Item Code :** 1220.2008.01

**Item Description :**

This Set Comprises Masses Of Brass With Slots, Brass Hangers With Hook. The Design Of The Slots Ensures That The Masses May Be Easily Removed And Replaced But Will Not Fall Of Accidentally.

- Set of 9 Slotted Masses, One Hanger Each Of 10gm – Total 100gm
- Set of 9 Slotted Masses, One Hanger Each Of 25gm – Total 250gm
- Set of 9 Slotted Masses, One Hanger Each Of 100gm – Total 1000gm
- Set of 4 Slotted Masses, One Hanger Each Of 50gm – Total 250gm
- Set of One Hanger 20gm, Three Slotted Masses 20g One Mass 10gm And Two Masses 5gm – Total 100gm
- Set of One Hanger 50gm, Nine Slotted Masses 20gm, One Mass Of 10gm And Two Masses Of 5gm – Total 250gm.



**Item Image :**

Cat. No.

1220.2008.01