

Item Name : Physical Balancee

Item Code : 3090.2008.01

Item Description :

It is sensitive precision balance.

Technical Data:

- Beam: Finest light brass beam
- Knife-edges and planes: made of finest agate.
- Arrestment: To beam and suspensions highly polished metal pillar fitted with six point Arrestment and extended support and single block guide.
- Case: wooden polished case with divided doors.
- Maximum Capacity: 250 Gms
- Sensitivity: 1mg.

Item Image :



Cat. No.

3090.2008.01