

Item Name : Precision balance

Item Code : 3101.2008.01

Item Description :

- Inbuilt motorized calibration (Optional)
- Multiple weighing units. • Standard bi-directional RS-232 interface
- Poly function : % weighing, counting, fill mode, GSM etc
- VFD display • Density determination Kit (optional)
- GSM (grams per sq mtr) measurement for paper & cloth quality checking.

Tech Specification:-

- Capacity :- 1000g • Readability :- 0.01g • Repeatability (+/-) :- 0.01g
- Linearity(+/-) :- 0.02g • Pan Size :- 125 mm

Item Image :



Cat. No.

3101.2008.01