

**Item Name :** Ripple tank apparatus

**Item Code :** 3109.2008.01

**Item Description :**

This economy Ripple Tank System comes as a ready-to-use unit with no external power source required. Comprised of a tough, moulded frame with a large transparent viewing surface. The Ripple Tank light source is a 6V halogen lamp which gives excellent illumination. The tank has sloping sides which minimise noise.

This Ripple Tank System comes complete with:

- Ripple motor, eccentric cam mounted on ripple bar
- Ripple tank
- 3 Support pillars for supporting the light source and rippler bar
- 4 Push-fit legs
- 2 Motor support springs
- Hand stroboscope
- 2 Ripplerbarhanging hooks
- Lamp unit with adjustable arm
- Concave perspex plate
- Rectangular perspex plate
- Curved reflector
- 2 Barriers
- Roller Bar

**Item Image :**



Cat. No.

3109.2008.01