

**Item Name :** Sodium Lamp

**Item Code :** 3112.2008.01

**Item Description :**

Low Pressure Sodium Vapor Lamp emits light at 589 nm and 589.6 nm, useful for wavelength calibration or resolution regulation. When the low-pressure Sodium lamp requires a warm up period of about ten minutes before it reaches substantially full light output, and has to cool down before it can be used again.

**Specifications**

- Power : 20 W • Operating Voltage : 20 VAC
- Operating Current : 1 ~ 1.3 A • Base Type : E27
- Dimension of Bulb : 28 mm in OD, 155 mm in length
- Dimension of Housing : 50 mm dia, 195 mm in length
- Dimension of Windows : 2 opposite windows: 30 mm height x 16 mm width
- Range of Height Adjustment : 140 mm to 240 mm

**Item Image :**



Cat. No.

3112.2008.01