

**Item Name :** Sonometer Wooden

**Item Code :** 3114.2008.01

**Item Description :**

It is a device developed in order to study vibrating strings.

Investigation on the dependency of the pitch on length, tension and thickness of a vibrating string can be done. Provided with three wires, two scales (scale length 50cm), one number of spring balance (capacity 100N), two fixed and three moveable bridges. A pulley is also attached on the corner of sonometer. Over all size 67.5x10x10.5cm approx.

**Item Image :**



Cat. No.

3114.2008.01