

Item Name : Voltmeter

Item Code : 3132.2008.01

Item Description :

This instruments for measuring voltage resp. current in shock-proof desktop housings built especially student and practical lab experiments with moving coil instruments, mirror scale and 4-mm safety sockets.

Specifications: • Measurement ranges: 3 V, 15 V, 300 V

• Scale division: 0.1 V, 1 V, 10 V • Internal resistance: 50 Ohm, 50 mV

• Accuracy: class 2.0 • Dimensions: 90x106x103 mm

Item Image :



Cat. No.

3132.2008.01