

Item Name : Resistance Box Plug Key

Item Code : 3133.2008.01

Item Description :

Teak wood box fitted with 7/16" rectangular brass, blocks accuracy +

0.1%. Provided resistance box is designed and manufactured by our deft professionals using premium quality components and advanced multiple resistors & switches to ensure appropriate flow of current, by adjusting the resistor value. This resistance box is available in the m


Features:

- Hassle-free operations • Long service life • Anti-corrosive body
- High durability • Sturdiness • Perfectly designed

Available In Range

- 1 to 10 or 1-1 • 1 to 50 or 1-5 • 1 to 100 or 1-10
- 1 to 500 or 1-50 • 1 to 1000 or 1-100 • 1 to 5000 or 1-500
- 1 to 10000 or 1-1000

Item Image :



Cat. No.

3133.2008.01