

**Item Name :** Semiconductor Diode Apparatus


**Item Code :** 3139.2008.01

**Item Description :**

Objective : To study Forward & Reverse Characteristics of PN Junction Diode.

Features : Instrument comprises of two DC Regulated Power Supplies 0-3VDC/ 150mA & 0-30 VDC/150mA, two round meters for voltage & current measurement, PN Junction Diodes mounted behind the panel, connections of Supplies, Meters & Diodes brought out at 4 mm Sockets.

**Item Image :**



Cat. No.

3139.2008.01