

Item Name : Compound Microscope

Item Code : 4009.2008.01

Item Description :

Body: Well designed heavy & sturdy die-casted monocular body inclinable to 90° suitable to meet the needs of college and school students. Tube with mechanical length 160mm. Focusing by separate coarse fine focusing knobs with triple revolving nose piece, horse shoe base. Square stage 120x120mm with clips for holding the slide. Double lens bright field ABBE Condenser is fixed with the stage with iris diaphragm and swing out filter holder. Illumination by plano concave mirror in fork mount.

Optical Combination: • Anti Fungal Objectives: 4/5x, 10x, 45x(SL)

• Eye Piece: 10x, 15x(Huygenian) • Packing: In thermocouple case with dust cover.

Item Image :



Cat. No.

4009.2008.01