

**Item Name :** Magnetising and Demagnetising Coil

**Item Code :** 3141.2008.01

**Item Description :**

The Magnetizing and Demagnetizing Coils low voltage and is suitable for magnetizing and demagnetizing ordinary magnets, iron bars, strips, etc. It comprises of a solenoid wound with insulated copper wire which is mounted on a base complete with switch and 4mm terminals.

- Solenoid 250mm long x 35mm internal diameter • Demagnetizing by a.c. only
- Operating voltage 12V a.c. or d.c. at 6A • Magnetizing by a.c. or d.c.

**Item Image :**



Cat. No.

3141.2008.01